

ANS & Arterial Assessment Device
Training Manual



nexus ANS+

PATIENTS SETTINGS MANUAL SUPPORT SIGN OUT

ARTERY ASSESSMENT

ARTERIAL RISK ASSESSMENT



ARTERIAL STIFFNESS

PERIPHERAL ARTERY DISEASE

Blood Pressure Pulse Volume Analysis		Ankle Brachial Index			
Left PWV	Right PWV	Left Ankle	Left Ankle PV	Right Ankle	Right Ankle PV
1142 cm/s	1026 cm/s	1.32 normal	615 normal	1.27 normal	322 normal
Peripheral An	38	Systolic (mmHg)	Diastolic (mmHg)	Systolic (mmHg)	Diastolic (mmHg)
CASP (mmHg)	76	132	82	127	80
		Left Ankle Blood Pressure	Right Ankle Blood Pressure	Brachial Blood Pressure	

AUTONOMIC NERVOUS SYSTEM ASSESSMENT

ANS RISK ASSESSMENT



SECRETORY FUNCTION

AUTONOMIC OSCILLATION

Autonomic Nervous System Response & Test Results	Heart Rate Variability	Cardiac Autonomic Index (CAI)
LF Amplitude (ms)	4807	SPS (ms) 1
LF Latency (ms)	960	Vl (ms) 1.37
LF Baseline (ms)	2.06	E (ms) 1.32
HF Amplitude (ms)	1049	CV (ms) 1.17
HF Latency (ms)	822	
HF Baseline (ms)	2.04	
HF Baseline (ms)	1074	

PHYSICIAN QUICK GUIDE MONITORING OF CHRx

SYMBOL LEGEND:

nexus ANS+

nexus ANS+

TABLE OF CONTENTS

Introduction.....	2
Test Procedure	4
Test Methods.....	5
Who Should be Tested	6
Contraindications.....	7
Results	8
Arterial Stiffness	10
Description.....	11
Key Measurements.....	12
Graphics & Ranges	13
Peripheral Artery Disease (PAD).....	14
Description.....	15
Key Measurements.....	16
Graphics & Ranges	17
Sudomotor Function.....	18
Description.....	19
Key Measurements.....	20
Graphics & Ranges	21
Autonomic Regulation.....	22
Description.....	23
Key Measurements.....	24
Graphics & Ranges	25
Cardiac Autonomic Neuropathy (CAN).....	26
Description.....	27
Key Measurements.....	28
Graphics & Ranges	29
Symptom Worksheet	30
Superbill.....	31

INTRODUCTION



Nexus ANS+ is a neurovascular analysis device that performs a series of fast, non-invasive tests on your nerves and arteries.

Nexus ANS+ assesses the Autonomic Nervous System and the arteries. The **Autonomic Nervous System (ANS)** is the part of the nervous system that regulates key involuntary functions of the body, including the activity of the heart, the smooth muscles (including the muscles of the intestinal tract), and the glands. The ANS has two divisions: the **Sympathetic Nervous System**, which accelerates the heart rate, constricts blood vessels, and raises blood pressure, and the **Parasympathetic Nervous System**, which slows the heart rate, increases intestinal and gland activity, and relaxes sphincter muscles.

Nexus ANS+ aids in the identification and early detection of ANS and arterial dysfunction related to disorders and disease complications. Autonomic Nerve and artery dysfunction or damage are the most frequent complications of **Chronic Metabolic Diseases (CMD)**, such as diabetes, kidney disease, hepatitis, thyroid failure, and aging.

Autonomic Neuropathy occurs when the nerves that control involuntary bodily functions are damaged. This may affect blood pressure, temperature control, digestion, bladder function and even sexual function. The nerve damage interferes with the messages sent between the brain and other organs and areas of the Autonomic Nervous System, such as the heart, blood vessels and sweat glands.

Arterial Disease is a general term that refers to a number of diseases that affect the blood vessels that distribute oxygenated blood throughout the body. Arterial Disease can hinder the body's ability to transport blood to affected body region(s). Depending on the type and severity of the disease, it can result in limb amputation and even death.

Although the symptoms may imitate each other, it is imperative to distinguish between Autonomic Neuropathy and Arterial Disease because the treatments are quite different. Nexus ANS+ gives healthcare providers the information necessary to make this distinction and improve patient treatment outcomes with early detection and accurate assessment.

Nexus ANS+ performs non-invasive tests to assess the patient's risk of autonomic nervous system & arterial dysfunction related to disorders and disease complications.




Arterial Risk Assessment

- Skin Microcirculatory Disorder
- Arterial Stiffness
- Peripheral Artery Disease

Skin Microcirculation, the circulation of blood to the smallest blood vessels in the skin, is altered in patients with hypertension, diabetes, and obstructive arterial disease. Furthermore, abnormalities of microcirculation are generally accepted as early changes associated with diabetes. Eventual manifestations of altered microcirculation, such as neuropathy, are related to the duration and severity of diabetes. Nexus ANS+ gathers information to evaluate risk of Skin Microcirculatory Disorder using the non-invasive Galvanic Skin Response (GSR) test.

Arterial Stiffness occurs as a consequence of biological aging and Arteriosclerosis, which occurs when arteries become thick and stiff, sometimes restricting blood flow to your organs and tissues. Arterial Stiffness is associated with an increased risk of cardiovascular events such as heart attack and stroke, the two leading causes of death in the developed world. Depending on the cause, Arterial Stiffness may be treated and prevented. Nexus ANS+ evaluates Arterial Stiffness using Pulse Wave Velocity (PWV), Peripheral augmentation index (pAIx), and Central aortic Systolic Pressure (CASP).

Peripheral Artery Disease (PAD) is a common manifestation of atherosclerotic vascular disease where the arteries in your legs or arms are narrowed or blocked. Its incidence increases with age and in the presence of known cardiovascular risk factors (e.g., smoking and diabetes). People with PAD are at an increased risk of heart attack, stroke, poor circulation and leg pain. Nexus ANS+ evaluates PAD using the Ankle Brachial Index (ABI). The ABI compares your systolic blood pressure measured at your ankle with your systolic blood pressure measured at your arm. A low ABI can be a strong indicator of PAD and risk of circulatory problems.



ANS Risk Assessment

- Sudomotor Function
- Autonomic Regulation
- Cardiac Autonomic Neuropathy

Sudomotor (or sweat motor) Function is related to the nerve fibers controlling the activity of the sweat glands (the post sympathetic cholinergic nerve fibers or C-fibers). Sudomotor dysfunction (sweat dysfunction) is an early indicator of small fiber neuropathy. Traditional neurophysiologic measurements of sudomotor function include invasive testing such as thermoregulatory sweat testing (TST), quantitative sudomotor axon reflex testing (QSART), silicone impressions. Nexus ANS+ uses Galvanic Skin Response (GSR) to test Sudomotor Function, which is a non-invasive test.

Autonomic Regulation is the body's ability to maintain homeostasis (stability and balance) during internal and external stimuli. Autonomic Regulation is always functioning and we are often unaware of the important tasks it is performing. When the nerves that control Autonomic Regulation are damaged, Autonomic Dysfunction can develop. Autonomic Dysfunction can be temporary or chronic. Diabetes and Parkinson's disease are two examples of chronic conditions that can lead to Autonomic Dysfunction. Nexus ANS+ tests Autonomic Regulation through a combination of Heart Rate Variability (HRV) Assessment and Cardiac Autonomic Reflex Tests (CARTs).

Cardiovascular Autonomic Neuropathy (CAN) is a common form of autonomic neuropathy, causing abnormalities in heart rate control and central/peripheral vascular dynamics. CAN has a strong link to diabetes and can contribute to the development of a variety of severe conditions resulting in a multitude of complications and higher mortality risk. Like the Autonomic Regulation Assessment, Nexus ANS+ evaluates the Sympathetic and Parasympathetic Nervous Systems using CARTs to identify CAN risk factors.

TEST PROCEDURE

Nexus ANS+ provides **fast, accurate, and clear results**. Easy to learn and operate, the tests can be administered by any staff member.

Testing with the Nexus ANS+ is performed in roughly 7 minutes and is completely non-invasive.

The recording includes a baseline phase where the patient is relaxed and a testing phase where the patient is asked to perform basic breathing exercises and one active postural change. Scores are calculated based on several recorded parameters and allow for the fast and intuitive interpretation of color coded results.



1. Set up the patient.

Input the patient's information into the Nexus ANS+ System.

2. Connect Nexus ANS+ to Patient.

With the patient laying supine:

- Place disposable electrodes on the soles of the feet
- Place blood pressure cuffs on the arm and legs
- Place pulse oximeter on the finger

3. Collect Baseline Data.

First the device will collect baseline information while the patient is relaxed.

- BP Pulse Volume Analysis
- Ankle Brachial Index (ABI)
- Galvanic Skin Response & Metaboreflex Analysis
- Heart Rate Variability (HRV) Analysis

4. Perform Cardiac Autonomic Reflex Tests (CARTs).

First the device will collect baseline information while the patient is relaxed.


- Valsalva Maneuver – patient bears down and attempts to exhale with mouth and nostrils closed
- Deep Breathing – patient performs a deep breathing exercise
- Change in Posture – patient stands up from a seated position

5. Test is Complete.

Collected data is placed through a series of algorithms that produce scores that indicate risk levels for a variety of ANS and Artery related diseases and conditions.



Nexus ANS+ performs **4 non-invasive tests in 7 minutes**. Results can aid in the early detection of Autonomic Nervous System and Arterial conditions.



Sudomotor Test

Used to assess risk of Small Fiber Neuroptathy & Skin Microcirculatory Disorder

The Nexus ANS+ uses the galvanic skin response method which is based on the induced sudomotor axon reflex stimulations of the cholinergic sympathetic system using:


- Direct Current at a constant voltage
- Metaboreflex following the inflation of the ankle cuffs

It provides quantitative evaluation of the sweat response measured on disposable electrodes.

Photoplethysmography (PTG)

Used to assess risk of Arterial Stiffness & Autonomic Dysregulation

Nexus ANS+ uses Photoplethysmography (PTG) analysis derived from the pulse oximeter. PTG is an optical measurement technique using red and infrared light to measure peripheral blood volume changes. PTG uses the time domain and spectral analysis of the waveform recording during 2 minutes of the test to retrieve information regarding artery (capillaries) stiffness, and Heart Rate Variability (HRV) Analysis.



Cardiac Autonomic Reflex Tests (CARTs)

Used to assess risk of Cardiac Autonomic Neuropathy

After a Baseline ANS Evaluation is completed by assessing the Heart Rate Variability (HRV) analysis, the Cardiac Autonomic Reflex Tests (CARTs) are performed. CARTs consist of 3 tests to assess the autonomic functions of the Sympathetic and Parasympathetic Systems.


The tests are as follows:

- Valsalva Maneuver
- Deep Breathing
- Change in Posture

Ankle Brachial Measurement (ABI)

Used to assess risk of Peripheral Artery Disease

The Ankle-Brachial Index (ABI) is the ratio of the blood pressure at the ankle to the blood pressure in the upper arm. Compared to the arm, lower blood pressure in the leg is an indication of blocked arteries due to arterial disease. The ABI is calculated by dividing the systolic blood pressure at the ankle by the systolic blood pressure in the arm.



WHO SHOULD BE TESTED

Over 45,000,000 Americans have symptoms or risk factors that would justify testing with the Nexus ANS+.

Symptom-Based

Any patient suffering any one of the following symptoms could be a candidate for the test:

- Pain, numbness, tingling, or burning in the feet
- Leg pain after effort, claudication
- Dizziness, Syncope

Age/Health-Based

Any patient above the age of 50 that has been diagnosed with any one of the following health conditions could be a candidate for the test:

- Hypertensive
- Smoker
- Overweight
- Diabetic

Age-Based

Any patient above the age of 70 could be a candidate for a test.



40% of adults age 50 and older are obese.



1 in 6 Americans aged 45-64 is a smoker.



Over 60% of Americans over the age of 55 have high blood pressure.



12 million Americans suffer from Lower Extremity Disease.



The CDC projects that 1 in 3 adults could have diabetes by 2050.



15-20% of people over 60 years old have Atery and ANS damage.

General Contraindications

The devices should not be used in association with or presence of defibrillators, cardiac pacemakers, patients connected to electronic life support devices, or any implanted electronic device.

Any blood pressure recording can be affected by the position of the subject, his or her physiologic condition, and other factors such as:

- Taking the blood pressure measurement:
 - » Less than 1 hour after meals
 - » After drinking alcohol or coffee
 - » Less than 1 hour after smoking
 - » After exercise or sporting activity
 - » After bathing
- White-Coat Hypertension (anxiety)
- Factors during Readings: Something as simple as moving the arm while blood pressure is being taken can affect the reading. If the sleeves of a shirt fit too tightly over the cuff, or if the cuff size is wrong, the readings could be inaccurate.
- Features: When the blood pressure reading is being taken, patient should be lying in a comfortable position.
- Considerations: Other factors include age, race, temperature, medication, emotional state, stress, exercise, smoking, and medical conditions such as diabetes, obesity, cardiovascular disease, irregular heartbeat, pregnancy and alcoholism. These factors need to be taken into account to determine an accurate measurement of blood pressure.
 - » The accuracy of the blood pressure device should be verified, and we recommend that a calibration is performed annually.
 - » The blood pressure device might not meet its performance specifications if stored or used outside of the temperature and humidity range.
 - » Blood pressure measurements determined with this device are equivalent to those obtained by an intra-arterial blood pressure measurement device within the limits prescribed by the American National Standard for manual, electronic, or automated sphygmomanometers.
 - » Cardiac arrhythmia may cause an irregular heartbeat, and may increase the measurement time.
 - » If the patient is on a Heart-Lung machine, the measurement may not be possible.
 - » Rapid pressure changes may not be possible to record.
 - » Severe shock or hypothermia may give unreliable results since reduced blood flow to the peripheries will reduce pulsation of arteries

Electrode Contraindications

- Dermatological lesions or calluses in contact with the electrodes, or excessive perspiration.
- Metal pins or prostheses on the level of the extremities or the joints.
- This device should not be used on pregnant women.
- An absence of one or more limbs.

ABI Contraindications

- Arterial catheter, AV fistula, or pressure dressing.
- Venous pulsations may cause erroneous readings in blood pressure (e.g. tricuspid valve regurgitation).
- Bilateral mastectomy

Oximeter Contraindications


- When using the oximeter probe, use the finger of the arm not in use with the blood pressure cuff, arterial catheter, or having an AV fistula or pressure dressing.
- Dyes introduced into the bloodstream, such as methylene blue, indocyanine green, indigo carmine, and fluorescein.
- Significant levels of dysfunctional hemoglobins, such as carboxyhemoglobin or methemoglobin, will affect the accuracy of the SpO2 measurement.
- Any condition that restricts blood flow, such as the use of a blood pressure cuff or extremes in systemic vascular resistance, may cause an inability to determine accurate pulse rate and SpO2 readings.
- For elderly patients or subjects with a weak pulse due to shock, low ambient/body temperature, major bleeding, or use of a vascular contracting drug, the SpO2 waveform will be decreased or absent.
- Fingernail polish or false fingernails may cause inaccurate SpO2 readings.
- Venous pulsations may cause erroneous readings in blood pressure (e.g. tricuspid valve regurgitation).
- Be careful with low perfused patients. Using the blood pressure device may cause skin erosion and/or pressure necrosis.
- Valsalva maneuver should not be performed on persons: Undergoing procedures with proliferative retinopathy.
- Systolic Blood pressure of 160mmHg or higher.
- Anyone who's had laser treatment for retinopathy within the past three months.

RESULTS



ARTERY ASSESSMENT

ARTERIAL RISK ASSESSMENT



ARTERIAL STIFFNESS

Blood Pressure Pulse Volume Analysis

Left PWV **1142** cm/s
Right PWV **1026** cm/s

Peripheral Alx **38**
CASP (mmHg) **76**

PERIPHERAL ARTERY DISEASE

Ankle Brachial Index

Left ABI **1.32** NORMAL
Left Ankle PV **619** mL/min
Right ABI **1.27** NORMAL
Right Ankle PV **322** mL/min

	Systolic (mmHg)	Diastolic (mmHg)	MAP (mmHg)
Left Ankle Blood Pressure	132	82	98
Right Ankle Blood Pressure	127	80	95
Brachial Blood Pressure	100	62	74

SYMBOL LEGEND:

Graphics/Ranges

Trend Chart

Off Label Use

- PATIENTS
- SETTINGS
- MANUAL
- SUPPORT
- SIGN OUT

AUTONOMIC NERVOUS SYSTEM ASSESSMENT

ANS RISK ASSESSMENT



SUDOMOTOR FUNCTION

Galvanic Skin Response & Metaboreflex

LF Amplitude (mV)	960
LF Latency (sec)	2.06
LF Baseline (mV)	1049
RF Amplitude (mV)	922
RF Latency (sec)	2.04
RF Baseline (mV)	1074

AUTONOMIC REGULATION

Heart Rate Variability

Total Power (ms ²)	4807
Stress Index (%)	38
HF (ms ²)	1782
SDANN (ms)	77

CAN


Cardiac Autonomic Reflex Tests


SPRS (mmHg)	1
Valsalva R.	1.37
E / I R.	1.32
K3015 R.	1.17

PHYSICIAN QUICK GUIDE

MONITORING OF CMR *

ARTERIAL STIFFNESS - RESULTS





Description

Arterial Stiffness occurs as a consequence of biological aging and Arteriosclerosis, which occurs when arteries become thick and stiff, sometimes restricting blood flow to your organs and tissues. Arterial Stiffness is associated with an increased risk of cardiovascular events such as heart attack and stroke, the two leading causes of death in the developed world. Depending on the cause, Arterial Stiffness may be treated and prevented. Nexus ANS+ evaluates Arterial Stiffness using Pulse Wave Velocity (PWV), Peripheral augmentation index (pAlx), and Central aortic Systolic Pressure (CASP).

Relevant Definitions

Arteriosclerosis occurs when the blood vessels that carry oxygen and nutrients from your heart to the rest of your body (arteries) become thick and stiff — sometimes restricting blood flow to your organs and tissues. Healthy arteries are flexible and elastic, but over time, the walls in your arteries can harden, a condition commonly called hardening of the arteries.

Atherosclerosis is a specific type of arteriosclerosis, but the terms are sometimes used interchangeably. Atherosclerosis refers to the buildup of fats, cholesterol and other substances in and on your artery walls (plaques), which can restrict blood flow. These plaques can burst, triggering a blood clot. Although atherosclerosis is often considered a heart problem, it can affect arteries anywhere in your body. Atherosclerosis may be preventable and is treatable.



Arterial Stiffness – Key Measurements

Pulse Wave Velocity

Pulse Wave Velocity (PWV) is a measure of arterial stiffness, or the rate at which pressure waves move down the vessel.

$$PWV = \frac{D \text{ (meters)}}{\Delta T \text{ (seconds)}}$$

A = Right Brachial Pulse pressure waveform

B = Right Ankle Pulse pressure waveform


ΔT = Delta Time between the 2 waves

D = Estimated distance from arm to ankle based on patient total height

Peripheral Augmentation Index

The Peripheral Augmentation Index (pAIx) is a ratio calculated from the blood pressure waveform, it is a measure of wave reflection and arterial stiffness. Peripheral Augmentation index is commonly accepted as a measure of the enhancement (augmentation) of central aortic pressure by a reflected pulse wave.

$$pAIx = \frac{P1}{P2} * 100$$



Pulse Pressure

Pulse Pressure (PP) is the difference between the systolic and diastolic pressure readings.

$$PP = SP - DP$$

SP = Systolic Pressure


DP = Diastolic Pressure

Central Aortic Systolic Pressure (CASP)

Central Aortic Systolic Pressure (CASP) is the blood pressure at the root of the aorta. It has been shown in many studies to be the strongest independent risk factor for stroke, heart attacks and the likelihood of survival after such an event.

$$CASP = DP + \frac{(pAIx * PP)}{100}$$

Arterial Stiffness – Graphics & Ranges



Measurement Ranges

PWV below 8.5m/s and/or Alx below 75%	PWV from 8.5 to 10m/s and/or Alx from 75 to 100%	PWV above 10m/s and/or Alx above 100 %
NORMAL: Suggests that the patient does not have arterial stiffness.	BORDERLINE: Suggests condition could be reversed by a lifestyle change.	ABNORMAL: Suggests moderate to severe arterial stiffness.

Clinical Data

Pressure measured with a cuff in the brachial artery is accepted as an important predictor of future cardiovascular risk. However, systolic pressure varies throughout the arterial tree, such that aortic (central) systolic pressure is actually lower than corresponding brachial values, although this difference is highly variable between individuals. Emerging evidence now suggests that central pressure is better related to future cardiovascular events than is brachial pressure. Moreover, anti-hypertensive drugs can exert differential effects on brachial and central pressure. Therefore, basing treatment decisions on central, rather than brachial pressure, is likely to have important implications for the future diagnosis and management of hypertension. Such a paradigm shift will, however, require further, direct evidence that selectively targeting central pressure, brings added benefit, over and above that already provided by brachial artery pressure. [McEnery C.M et al.2014]


Educational Reading

Determination of age-related increases in large artery stiffness by digital pulse contour analysis.
 Department of Clinical Pharmacology, St. Thomas' Hospital, Centre for Cardiovascular Biology and Medicine, King's College London, Lambeth Palace Road, London SE1 7EH, UK.
<http://www.ncbi.nlm.nih.gov/pubmed/12241535>

M. Span¹, G. Gersak², M. Meza³, A. Kosir⁴

https://asiakas.kotisivukone.com/files/terveystekniikka.palvelee.fi/Tutkimukset_artikkelit/mesi_clinical_study.pdf

PERIPHERAL ARTERY DISEASE (PAD) - RESULTS



Description

Peripheral Artery Disease (PAD) is a common manifestation of atherosclerotic vascular disease where the arteries in your legs or arms are narrowed or blocked. Its incidence increases with age and in the presence of known cardiovascular risk factors (e.g., smoking and diabetes). People with PAD are at an increased risk of heart attack, stroke, poor circulation and leg pain. Nexus ANS+ evaluates PAD using the Ankle Brachial Index (ABI). The ABI compares your systolic blood pressure measured at your ankle with your systolic blood pressure measured at your arm. A low ABI can be a strong indicator of PAD and risk of circulatory problems.

PERIPHERAL ARTERY DISEASE			
Ankle Brachial Index			
Left ABI	Left Ankle PV	Right ABI	Right Ankle PV
1.32 NORMAL	619 mL/min	1.27 NORMAL	322 mL/min
	Systolic (mmHg)	Diastolic (mmHg)	MAP (mmHg)
Left Ankle Blood Pressure	132	82	98
Right Ankle Blood Pressure	127	80	95
Brachial Blood Pressure	100	62	74

Peripheral Artery Disease (PAD) – Key Measurements

Ankle Brachial Index (ABI) is the ratio of the systolic pressure of the ankle divided by the arm systolic pressure.

Pulse Volume Velocity (PVV) reflects the return of the blood volume over time after inflation and occlusion.

Mean Arterial Pressure (MAP) is the average pressure in the arteries during one cardiac cycle.


Ankle Brachial Index (ABI)

The ankle-brachial index test is a quick, noninvasive way to check your risk of peripheral artery disease (PAD). Peripheral artery disease is a condition in which the arteries in your legs or arms are narrowed or blocked. People with peripheral artery disease are at an increased risk of heart attack, stroke, poor circulation and leg pain.

The ankle-brachial index test compares your systolic blood pressure measured at your ankle with your systolic blood pressure measured at your arm. A low ankle-brachial index number can indicate narrowing or blockage of the arteries in your legs, increasing your risk of circulatory problems, and possibly causing heart disease or stroke.

The Ankle Brachial Index (ABI) was initially proposed for the noninvasive diagnosis of lower-extremity peripheral artery disease (PAD) (Carter 1967). However, it was later shown that the ABI is also an indicator of atherosclerosis at other vascular sites, and it can serve as a prognostic marker for cardiovascular events and functional impairment, even in the absence of symptoms of PAD (Criqui et al. 1992).

The ABI results interpretation come from the “Guidelines for the Management of Patients with Peripheral Artery Disease” published by the Journal of the American College of Cardiology.



Peripheral Artery Disease (PAD) - Graphics & Ranges

PVR GRAPHICS RECORDS ✕

PULSE VOLUME RECORDS

Left Ankle PVR

Right Ankle PVR

Arm PVR

MARKERS RANGE SCALE

ABI

Severe	Abnormal	Borderline	Normal	NC
0.6	0.9	1.0	1.4	

Pulse Volume

Abnormal	Borderline	Normal
80	120	

An ABI less than 0.90 has been shown to have a sensitivity of 90% and a specificity of 98% for detecting a lower-extremity stenosis of greater than 50%. (Park C. W. 2013)

Measurement Ranges

ABPI less than 0.6	ABPI from 0.6 to 0.9	ABPI from 0.9 to 0.99	ABPI from 1.0 to 1.3	ABPI from 1.31 to 1.39	ABPI more than 1.4
SEVERE: Blockage in ankle and leg arteries may be severe.	ABNORMAL: Blockage in ankle and leg arteries may be fairly significant.	BORDERLINE: Ankle and leg arteries may be narrowing.	NORMAL: Ankle and leg arteries are likely normal.	BORDERLINE: Ankle and leg arteries may be borderline rigid.	ABNORMAL: Ankle and leg arteries may be rigid and do not compress.

PVR Max less than 50 mL	PVR Max from 50 to 79 mL	PVR Max from 80 to 119 mL	PVR Max more than 120 mL
Severe reduction in blood flow of the foot.	Moderate reduction in blood flow of the foot.	Borderline reduction in blood flow of the foot.	Normal blood flow of the foot.

Clinical Data

When comparing ABI measurement results obtained by oscillometric measurements vs. Doppler probe measurements, it was evident that the oscillometric measurements were more precise than Doppler probe measurements; 17% vs. 19% respectively. Bland Altman plot showed bias of +0.06 and scatter diagram showed good correlation (slope: 0.75). In the oscillometric method analysis percentage error was 0.26. Furthermore, assessment of the concordance index revealed a very good agreement between both methods in terms of clinical relevance: concordance index of 0.88 (95% CI: from 0.76 to 0.97). Oscillometric measurements were completed three times faster than Doppler probe measurements.


Educational Reading

Comparison of Doppler probe method and oscillometric method in determining ankle brachial index

M. Span¹, G. Gersak², M. Meza³, A. Kosir⁴

https://asiakas.kotisivukone.com/files/terveystekniikka.palvelee.fi/Tutkimukset_artikkelit/mesi_clinical_study.pdf

SUDOMOTOR FUNCTION - RESULTS



SUDOMOTOR FUNCTION	
Galvanic Skin Response & Metaboreflex	
LF Amplitude (mV)	960
LF Latency (sec)	2.06
LF Baseline (mV)	1049
RF Amplitude (mV)	922
RF Latency (sec)	2.04
RF Baseline (mV)	1074

Description

Sudomotor (or sweat motor) Function is related to the nerve fibers controlling the activity of the sweat glands (the post sympathetic cholinergic nerve fibers or C-fibers). Sudomotor dysfunction (sweat dysfunction) is an early indicator of small fiber neuropathy. Traditional neurophysiologic measurements of sudomotor function include invasive testing such as thermoregulatory sweat testing (TST), quantitative sudomotor axon reflex testing (QSART), silicone impressions. Nexus ANS+ uses Galvanic Skin Response (GSR) to test Sudomotor Function, which is a non-invasive test.

Relevant Definitions

Sudomotor (from Latin sudor, 'sweat' and motor) describes anything that stimulates the sweat glands.

Peripheral Neuropathy (PN) is damage to or disease affecting nerves, which may impair sensation, movement, gland or organ function, or other aspects of health, depending on the type of nerve affected.

Small Fiber Peripheral Neuropathy is a type of peripheral neuropathy that occurs from damage to the small unmyelinated peripheral nerve fibers. These fibers, categorized as C fibers, are present in skin, peripheral nerves, and organs. The role of these nerves is to innervate the skin (somatic fibers) and help control autonomic function (autonomic fibers). It is estimated that 15-20 million people in the United States suffer from some form of peripheral neuropathy.

Afferent nerves from skeletal muscle are stimulated by contraction, leading to increased ventilation and increased efferent sympathetic nerve activity. This reflex arc is called the muscle chemoreflex or **Metaboreflex**.

Galvanic Skin Response (GSR) is the phenomenon that the skin momentarily becomes a better conductor of electricity when either external or internal stimuli occur that are physiologically arousing.

Sudomotor Function - Key Measurements

Left/Right Foot Amplitude (L/RF Amplitude) reflects the sudomotor response to electrical and muscle metaboreflex stimulation of the sympathetic system. The response is proportionate to the active sweat gland density and skin microcirculation response. This is the main marker of sudomotor function.

Left/Right Foot Latency (L/RF Latency) reflects the time of sudomotor response after stimulation related to the skin thermoreceptor (A Delta Fiber).

Left/Right Foot Baseline (L/RF Baseline) is a marker of the left foot skin microcirculation.

Galvanic Skin Response (GSR)


The Nexus ANS + performs **galvanic skin response** based on the sympathetic skin response method and induced sudomotor axon reflex of the cholinergic sympathetic system using double stimulation:

- 1) DC electricity at a constant voltage
- 2) Metaboreflex following the inflation of the ankle pressure cuffs.


It provides a quantitative evaluation of the sweat response or sudomotor function measured via disposable electrodes. The quantitative induced sweat response is expressed by waveform amplitude from the skin voltage responses (in mV) and latency (in seconds).

The skin voltage responses reflect the absorption of the induced sweat on the bulk of the cloth electrodes, and latency reflects the time of response from stimulation to the peak of waveform amplitude. The measurement is performed at the soles of the right and left foot, which have the highest density of eccrine sweat glands and C-fibers that travel the furthest from the brain and spinal cord.

As perspiration increases, more sweat glands are stimulated which increase the voltage waveform amplitude in a given area of skin covered by the disposable cloth electrodes.



Sudomotor Function - Graphics & Ranges



Clinical Data


Sudomotor response (SMR) score. SMR score is based on highest score between the right and left foot. SMR Score had a sensitivity of 91.4% and specificity of 79.1% to detect peripheral distal neuropathy symptoms (Gandhi et al. 2015). The sudomotor testing clinical data suggests it may be the most sensitive means to detect peripheral small fiber neuropathy (Low et al. 2006). Amplitude response reflects the C-Fiber density, Baseline response reflects the skin microcirculation state and latency response reflects the small myelinated thermal fiber (A delta fiber) state following the induced sudomotor axon reflex stimulations of the cholinergic sympathetic system by electrical stimulation and metaboreflex.


Educational Reading

Detection of neuropathy using a sudomotor test in type 2 diabetes
Gandhi PG, Rao GHR
https://www.dovepress.com/articles.php?article_id=19898

Cardiac autonomic neuropathy in patients with diabetes mellitus
Gerasimos Dimitropoulos, Abd A Tahrani, and Martin J Stevens
<http://www.ncbi.nlm.nih.gov/pmc/articles/PMC3932425/>

RESULTS - AUTONOMIC REGULATION





Description

Autonomic Regulation is the body's ability to maintain homeostasis (stability and balance) during internal and external stimuli. Autonomic Regulation is always functioning and we are often unaware of the important tasks it is performing. When the nerves that control Autonomic Regulation are damaged, Autonomic Dysfunction can develop. Autonomic Dysfunction can be temporary or chronic. Diabetes and Parkinson's disease are two examples of chronic conditions that can lead to Autonomic Dysfunction. Nexus ANS+ tests Autonomic Regulation through a combination of Heart Rate Variability (HRV) Assessment and Cardiac Autonomic Reflex Tests (CARTs).

Autonomic Regulation – Key Measurement Details

Total Power is the main indicator of ANS activity at rest. In peer reviews, Total Power is used to define the total HRV Score.

High Frequency (HF) of the RR Intervals spectral analysis. It reflects Parasympathetic activities at rest.

Stress Index reflects the sympathetic system's effect on peripheral circulation. Sympathetic effect of peripheral small fiber.

Standard Deviation Average of all Normal to Normal (SDANN) is an indicator of both sympathetic and parasympathetic regulation at rest, and is a marker of V02 max.

Heart Rate Variability (HRV)

Heart Rate Variability (HRV) is the variation in the time interval between heartbeats. It is measured by the variation in the beat-to-beat intervals. It's measured in milliseconds. The Autonomic Nervous System plays a central role in HRV. Depressed or reduced HRV primarily means a lowered ability of the ANS regulatory function to keep homeostasis, cope with internal and external stressors (stress provoking agents), resist disease and/or recover in proper time.


Photoplethysmography (PTG)

Photoplethysmography (PTG) is an optical measurement technique using red and infrared light to measure peripheral blood volume obtained with the pulse oximeter.

- Is reflective of blood movement in cutaneous vessels and synchronous cardiovascular events.
- The signal contains information regarding the heart rate and variability, vessel dilation and contraction.



Autonomic Regulation – Graphics & Ranges



Clinical Data

A reduction of HRV and High score of Autonomic Regulation Score have been associated with the early stages of Cardiac Autonomic Neuropathy (CAN). HRV analysis reflects the autonomic nervous system components balance and activity. HRV results could be improved by lifestyle change.

Educational Reading

Heart rate variability: Standards of measurement, physiological interpretation, and clinical use: Task Force of The European Society of Cardiology and The North American Society of Pacing and Electrophysiology (Membership of the Task Force listed in the Appendix) https://www.escardio.org/static_file/Escardio/Guidelines/Scientific-Statements/guidelines-Heart-Rate-Variability-FT-1996.pdf


Impact of reduced heart rate variability on risk for cardiac events. The Framingham Heart Study.

Hisako Tsuji, MD; Martin G. Larson, ScD; Ferdinand J. Venditti, Jr, MD; Emily S. Manders, BS; Jane C. Evans, MPH; Charles L. Feldman, ScD; Daniel Levy, MD
<http://www.ncbi.nlm.nih.gov/pubmed/8941112>

Photoplethysmography and its application in clinical physiological measurement.

Regional Medical Physics Department, Freeman Hospital, Newcastle upon Tyne, UK. john.allen@nuth.nhs.uk
<http://www.ncbi.nlm.nih.gov/pubmed/17322588>

RESULTS - CARDIAC AUTONOMIC NEUROPATHY (CAN)



CAN	
Cardiac Autonomic Reflex Tests	
SPRS (mmHg)	1
Valsalva R.	1.37
E / IR.	1.32
K3015 R.	1.17

Description

Cardiovascular Autonomic Neuropathy (CAN) is a common form of autonomic neuropathy, causing abnormalities in heart rate control and central/peripheral vascular dynamics. CAN has a strong link to diabetes and can contribute to the development of a variety of severe conditions resulting in a multitude of complications and higher mortality risk. Like the Autonomic Regulation Assessment, Nexus ANS+ evaluates the Sympathetic and Parasympathetic Nervous Systems using CARTs to identify CAN risk factors.

Autonomic Regulation (CAN) – Key Measurements

Systolic Pressure Response to Standing (SPRV) is the difference between the sitting systolic pressure minus the standing systolic pressure. It reflects the Sympathetic adrenergic response.

Valsalva Ratio (VALS. R) is the subsequent tachycardia/bradycardia during the Valsalva maneuver are baro receptor mediated.

Expiration and Inspiration Ratio (E/I R) is calculated by the longest RR or NN Interval divided by the shortest RR Interval during the deep breathing test. It reflects the parasympathetic system's cardio-vagal response during breathing challenge.

K3015 R is the maximum RR intervals at 30 seconds divided by the minimum RR intervals at 15 seconds after standing up. It reflects the parasympathetic system cardiovagal response to the changes of posture. It is considered an index of cardiovascular function.

Cardiac Autonomic Reflex Tests (CARTs)

After a Baseline ANS Evaluation is completed by assessing the Heart Rate Variability (HRV) analysis, the **Cardiac Autonomic Reflex Tests (CARTs)** are performed. CARTs consist of 3 tests to assess the autonomic functions of the Sympathetic and Parasympathetic Systems. The tests are as follows:

Valsalva Maneuver – patient bears down and attempts to exhale with mouth and nostrils closed

Deep Breathing – patient performs a deep breathing exercise

Change in Posture – patient stands up from a seated position


Cardiac Autonomic Reflex Tests (CARTs) are recognized in order to diagnose Cardiac Autonomic Neuropathy (CAN) by the guidance of the Cardiac Autonomic Neuropathy (CAN) Subcommittee of the Toronto Consensus Panel on Diabetic Neuropathy and American Academy of Neurology. Cardiac Autonomic neuropathy occurs when the nerves that control involuntary bodily functions are damaged. This may affect blood pressure, temperature control, digestion, bladder function and even sexual function.

Cardiovascular (Parasympathetic) Evaluation with Valsalva ratio: this is the ratio of heart rate phase IV/heart rate phase II, and is essentially a measure of vagal function. The normal response is an increase in heart rate during phase II (SPVR2) in response to the fall in blood pressure, and the baroreflex response to the blood pressure overshoot in phase IV (SPRV4) is transient bradycardia. In dysautonomia patients, there is typically a loss of both the blood pressure overshoot and the reflex bradycardia. Age-adjusted normal Valsalva ratios are used as controls for the Valsalva ratio.

Adrenergic (Sympathetic) Evaluation with Beat-to-beat blood pressure response to the Valsalva maneuver: phase II may be divided into early and late phases. In early phase II, reduced preload and stroke volume lead to a fall in cardiac output despite the tachycardia caused by decreased vagal activity. Total peripheral resistance increases due to increased sympathetic discharge (nucleus tractus solitarius) and an increase in plasma epinephrine concentration. The fall in blood pressure is thus arrested, and is termed "late phase II"; in normal patients there is actually a rise in blood pressure just before release. Interestingly, this phase can be blocked with alpha-antagonists (e.g., phentolamine) while phase IV can be blocked with beta-blockers. Patients with peripheral adrenergic failure due to involvement of autonomic fibers will have absent late phase II and in fact can have such dramatic drops in blood pressure during forced expiration that syncope may occur.

Deep Breathing tests cardiac parasympathetic functions. Because the heart responses to deep breathing are mediated by the vagal nerve, the test is also referred to as cardiovascular testing. During inspiration, intra thoracic pressure is reduced, which will decrease venous return from the pulmonary circulation to the left atrium of the heart and decrease preload; heart rate will therefore increase slightly to compensate and maintain the same cardiac output. The opposite happens during expiration; intra thoracic pressure rises, increasing pulmonary venous return (and therefore preload), so heart rate will decrease to maintain constant cardiac output. The results is expressed as E/I Ratio (expiration/Inspiration Ratio).

Cardiac Autonomic Neuropathy (CAN) – Graphics & Ranges



Educational Reading

1. Cardiovascular Autonomic Neuropathy in Diabetes: clinical impact, assessment, diagnosis, and management
Spallone v, Ziegler o, Freeman R, Bernardi L, Frontoni s, Pop-Busul R, Stevens M, Kempler P, Hilsted J, Testaye s, Low P, vatenst P;
Toronto Consensus Panel on Diabetic Neuropathy. Diabetes Metab Res Rev. 2011 Oc~27(7):639-53. doi: 10.1002/dmrr.1239.

Diagnosis and management of diabetic autonomic neuropathy
"D J EWING, B F CLARKE 1982 BRITISH MEDICAL JOURNAL VOLUME 285"

PATIENT INFORMATION

Name: _____ **Date of Birth:** _____ **Gender:** Male Female
Height: _____ **Weight:** _____ **Activity Level:** None Light Moderate Frequent Heavy

Have you experienced any of the following in the last 14 days?

FOR OFFICE USE ONLY

			ANS	Vascular
Frequent Urination	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Nausea	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Numbness in Hands & Feet	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Vomiting	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Burning Sensation	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Difficulty Digesting Food	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Sweat Abnormalities	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Sexual Difficulties	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Tingling in Hands & Feet	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Urinary Problems	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Cold, Clammy, Pale Skin	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Depression	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Fainting	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Lack of Concentration	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Lack of Energy	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Rapid Shallow Breathing	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Exercise Intolerance	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Painful Contact with Sock or Bed	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Pebble or Sand-like Sensation in Shoes	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Stabbing or Electrical Shock Sensation	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Pins & Needles in Feet	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Blurred Vision	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Elevated Blood Sugar	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Extreme Thirst	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Fatigue	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Increased Hunger	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Heartburn	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Pain In Calves	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Angina	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Chest Pain that Goes Away with Rest	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Shortness of Breath	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
TIA (Mini Stroke)	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Blood Clot in a Vein	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Heart Attack	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Irregular Heartbeat (Too Fast/Slow)	<input type="checkbox"/> Yes	<input type="checkbox"/> No		
Stroke	<input type="checkbox"/> Yes	<input type="checkbox"/> No		

Patient Name: _____

Gender: Male Female

Date of Birth: _____

CODING INFORMATION AND SUGGESTED BILLING AMOUNTS

CODES & MODIFIERS

95923 - \$325.00 **95921** - \$150.00
93922 - \$180.00 **Modifier:** XU or 59

- CPT codes 95921 and 95943 CANNOT be used together.
- Use the 76 modifier if the physician has the need to repeat the test.
- Physicians can bill an E/M code along with these services. List E/M with a 25 Modifier then list test performed.
- Bill 1 unit for each CPT code
- Modifiers XU or 59 must be used on secondary CPT codes. Medicare allows modifier XU and 59 for Commercial payers.

CODES SPECIFIC TO 95921 AND 95923

<input type="checkbox"/> E08.40 Diabetes mellitus due to underlying condition with diabetic neuropathy, unspecified	<input type="checkbox"/> E85.0 Non-neuropathic hereditary amyloidosis
<input type="checkbox"/> E08.41 Diabetes mellitus due to underlying condition with diabetic mononeuropathy	<input type="checkbox"/> E85.1 Neuropathic hereditary amyloidosis
<input type="checkbox"/> E08.42 Diabetes mellitus due to underlying condition with diabetic polyneuropathy	<input type="checkbox"/> E85.2 Hereditary amyloidosis, unspecified
<input type="checkbox"/> E08.43 Diabetes mellitus due to underlying condition with diabetic autonomic	<input type="checkbox"/> E85.3 Secondary systemic amyloidosis
<input type="checkbox"/> E08.44 Diabetes mellitus due to underlying condition with diabetic amyotrophy	<input type="checkbox"/> E85.4 Organ-limited amyloidosis
<input type="checkbox"/> E08.49 Diabetes mellitus due to underlying condition with other diabetic neurological complication	<input type="checkbox"/> E85.8 Other amyloidosis
<input type="checkbox"/> E08.610 Diabetes mellitus due to underlying condition with diabetic neuropathic arthropathy	<input type="checkbox"/> E85.9 Amyloidosis, unspecified
<input type="checkbox"/> E08.65 Diabetes mellitus due to underlying condition with hyperglycemia	<input type="checkbox"/> G90.1 Disorders of autonomic nervous system, Familial dysautonomia
<input type="checkbox"/> E10.40 Type 1 diabetes mellitus with diabetic neuropathy, unspecified	<input type="checkbox"/> G90.2 Horner's syndrome
<input type="checkbox"/> E10.41 Type 1 diabetes mellitus with diabetic mononeuropathy	<input type="checkbox"/> G90.3 Multi-system degeneration of the autonomic nervous system
<input type="checkbox"/> E10.42 Type 1 diabetes mellitus with diabetic polyneuropathy	<input type="checkbox"/> G90.4 Autonomic dysreflexia
<input type="checkbox"/> E10.43 Type 1 diabetes mellitus with diabetic autonomic (poly)neuropathy	<input type="checkbox"/> G90.50 Complex regional pain syndrome I, unspecified
<input type="checkbox"/> E10.44 Type 1 diabetes mellitus with diabetic amyotrophy	<input type="checkbox"/> G90.511 Complex regional pain syndrome I of right upper limb
<input type="checkbox"/> E10.49 Type 1 diabetes mellitus with other diabetic neurological complication	<input type="checkbox"/> G90.512 Complex regional pain syndrome I of left upper limb
<input type="checkbox"/> E10.610 Type 1 diabetes mellitus with diabetic neuropathic arthropathy	<input type="checkbox"/> G90.513 Complex regional pain syndrome I of upper limb, bilateral
<input type="checkbox"/> E10.65 Type 1 diabetes mellitus with hyperglycemia	<input type="checkbox"/> G90.519 Complex regional pain syndrome I of unspecified upper limb
<input type="checkbox"/> E13.40 Other specified diabetes mellitus with diabetic neuropathy, unspecified	<input type="checkbox"/> G90.521 Complex regional pain syndrome I of right lower limb
<input type="checkbox"/> E13.41 Other specified diabetes mellitus with diabetic mononeuropathy	<input type="checkbox"/> G90.522 Complex regional pain syndrome I of left lower limb
<input type="checkbox"/> E13.42 Other specified diabetes mellitus with diabetic polyneuropathy	<input type="checkbox"/> G90.523 Complex regional pain syndrome I of lower limb, bilateral
<input type="checkbox"/> E13.43 Other specified diabetes mellitus with diabetic autonomic (poly)neuropathy	<input type="checkbox"/> G90.529 Complex regional pain syndrome I of unspecified lower limb
<input type="checkbox"/> E13.44 Other specified diabetes mellitus with diabetic amyotrophy	<input type="checkbox"/> G90.59 Complex regional pain syndrome I of other specified site
<input type="checkbox"/> E13.49 Other specified diabetes mellitus with other diabetic neurological complication	<input type="checkbox"/> G90.8 Other disorders of autonomic nervous system
<input type="checkbox"/> E11.40 Type 2 diabetes mellitus with diabetic neuropathy, unspecified	<input type="checkbox"/> G90.9 Disorder of the autonomic nervous system, unspecified
<input type="checkbox"/> E11.41 Type 2 diabetes mellitus with diabetic mononeuropathy	<input type="checkbox"/> G60.0 Hereditary motor and sensory neuropathy
<input type="checkbox"/> E11.42 Type 2 diabetes mellitus with diabetic polyneuropathy	<input type="checkbox"/> G60.1 Refsum's disease
<input type="checkbox"/> E11.43 Type 2 diabetes mellitus with diabetic autonomic (poly)neuropathy	<input type="checkbox"/> G60.2 Neuropathy in association with hereditary ataxia
<input type="checkbox"/> E11.44 Type 2 diabetes mellitus with diabetic amyotrophy	<input type="checkbox"/> G60.3 Idiopathic progressive neuropathy
<input type="checkbox"/> E11.49 Type 2 diabetes mellitus with other diabetic neurological complication	<input type="checkbox"/> G60.8 Other hereditary and idiopathic neuropathies
<input type="checkbox"/> E11.610 Type 2 diabetes mellitus with diabetic neuropathic arthropathy	<input type="checkbox"/> G60.9 Hereditary and idiopathic neuropathy, unspecified
<input type="checkbox"/> E11.65 Type 2 diabetes mellitus with hyperglycemia	<input type="checkbox"/> G61.81 Chronic inflammatory demyelinating polyneuritis
<input type="checkbox"/> E13.610 Other specified diabetes mellitus with diabetic neuropathic arthropathy	<input type="checkbox"/> G63 Polyneuropathy in diseases classified elsewhere
<input type="checkbox"/> R55 Syncope and collapse	<input type="checkbox"/> I49.8 Other specified cardiac arrhythmias
<input type="checkbox"/> R00.0 Tachycardia, unspecified	

CODES SPECIFIC TO 93922

<input type="checkbox"/> I72.4 Aneurysm Of Artery Of Lower Extremity.	<input type="checkbox"/> E10.65 Diabetes Mellitus W/ P Circ Disorders I Uncontrolled
<input type="checkbox"/> I72.1 Aneurysm Of Artery Iliac Crest	<input type="checkbox"/> E11.51 Diabetes Mellitus W/ P Circ Disorders II Not Stated As Uncontrolled.
<input type="checkbox"/> I72.3 Aneurysm Of Artery Upper Extremity.	<input type="checkbox"/> E11.65 Diabetes Mellitus W/ P Circ Disorders, II Unspecif Type Uncontrolled.
<input type="checkbox"/> I70.209 Atherosclerosis Native Arteries Extremities Unspecified.	<input type="checkbox"/> E08.51 Diabetes mellitus due to underlying condition W/ diabetic peripheral angiopathy W/O gangrene
<input type="checkbox"/> I70.229 Atherosclerosis Native Arteries Extremities With Rest Pain.	<input type="checkbox"/> E08.52 Diabetes mellitus due to underlying condition W/ diabetic peripheral angiopathy W/ gangrene
<input type="checkbox"/> I70.269 Atherosclerosis Native Arteries W/ Gangrene	<input type="checkbox"/> E09.51 Drug or chemical induced diabetes mellitus W/ diabetic peripheral angiopathy W/O gangrene
<input type="checkbox"/> I70.219 Atherosclerosis Native Arteries W/ Intermittent Claudication.	<input type="checkbox"/> E09.52 Drug or chemical induced diabetes mellitus W/ diabetic peripheral angiopathy W/ gangrene
<input type="checkbox"/> I70.25 Atherosclerosis Native Arteries W/Ulceration.	<input type="checkbox"/> E10.51 Type 1 diabetes mellitus W/ diabetic peripheral angiopathy W/O gangrene
<input type="checkbox"/> I70.599 Atherosclerosis Non Autologous Biological Bypass Graft Extremities.	<input type="checkbox"/> E10.52 Type 1 diabetes mellitus W/ diabetic peripheral angiopathy W/ gangrene
<input type="checkbox"/> I70.0 Atherosclerosis Of Aorta	<input type="checkbox"/> E11.52 Type 2 diabetes mellitus W/ diabetic peripheral angiopathy W/ gangrene
<input type="checkbox"/> I70.499 Atherosclerosis Of Autologous Vein Bypass Graft Extremities.	<input type="checkbox"/> E13.52 Other specified diabetes mellitus W/ diabetic peripheral angiopathy W/ gangrene
<input type="checkbox"/> I70.399 Atherosclerosis Of Bypass Graft Extremities	<input type="checkbox"/> I10 Hypertension Unspecified
<input type="checkbox"/> I77.71 Dissection of carotid artery	<input type="checkbox"/> I30.9 Acute pericarditis, unspecified
<input type="checkbox"/> I77.72 Dissection of iliac artery	<input type="checkbox"/> I40.0 Infective myocarditis
<input type="checkbox"/> I77.73 Dissection of renal artery	<input type="checkbox"/> I40.1 Isolated myocarditis
<input type="checkbox"/> I77.74 Dissection of vertebral artery	<input type="checkbox"/> I48.0 Paroxysmal atrial fibrillation
<input type="checkbox"/> I67.0 Dissection of cerebral arteries, nonruptured	<input type="checkbox"/> I48.1 Persistent atrial fibrillation
<input type="checkbox"/> I96 Gangrene, not elsewhere classified	<input type="checkbox"/> I48.2 Chronic atrial fibrillation
<input type="checkbox"/> I73.9 Pheripheral Vascular Disease, Unspecified	<input type="checkbox"/> I48.91 Unspecified atrial fibrillation
<input type="checkbox"/> I73.00 Raynaud's syndrome without gangrene	<input type="checkbox"/> Z13.9 Screening For Hypertension
<input type="checkbox"/> I73.01 Raynaud's syndrome with gangrene	<input type="checkbox"/> G54.0 Brachial plexus disorders
<input type="checkbox"/> I73.1 Thromboangiitis obliterans (Buerger's disease)	<input type="checkbox"/> I74.01 Saddle embolus of abdominal aorta

Ordering Physician: _____ Signature: _____ Date: _____

* The information contained in this template superbill is intended as general information only. It is not intended to serve as a substitute for professional advice of a physician or medical coding professional. Final determination should be made by the treating physician. CODES SPECIFIC TO 93922, 95921 & 95923. * Physicians can bill an E/M code along with these services.



nexus ANS+

866.732.5119 • www.nexusans.com

Distributed by:  BYNIX